



MANIPAL INSTITUTE OF TECHNOLOGY
MANIPAL
(A constituent unit of MAHE, Manipal)

**DEPARTMENT OF DATA SCIENCE & COMPUTER
APPLICATIONS**

Medical Image Processing

Submitted by

Suzen Malik Firasta

210968058

Under the guidance of

Mr. Girisha Surathkal [MAHE-MIT]

Assistant Professor

Dept. of Data Science and Computer Applications

Advanced Medical Image Processing Tool

Suzen Malik Firasta

November 8, 2024

1 Introduction

This report outlines the functionalities and results of the Advanced Medical Image Processing Tool, implemented using Streamlit and OpenCV. The tool provides various image processing techniques for medical images, including quality measures, enhancements, noise reduction, edge detection, and segmentation.

Welcome to the **Advanced Medical Image Processing Tool**—an intuitive, powerful, and user-friendly application crafted to streamline medical image analysis for healthcare professionals, researchers, and enthusiasts alike. Using this tool, you can easily analyze and enhance medical images with just a few clicks, gaining deep insights and clarity for your images. Whether you're working with X-rays, MRIs, or any grayscale medical imaging, this tool is designed to make complex processes accessible and interactive.

2 Features

2.1 Quality Measures

- **Michelson Contrast:** Measures the contrast between the maximum and minimum pixel intensity in the image.

Formula:

$$C = \frac{I_{\max} - I_{\min}}{I_{\max} + I_{\min}} \quad (1)$$

- **RMS Contrast:** Measures the standard deviation of pixel intensities.

Formula:

$$C_{\text{RMS}} = \sqrt{\frac{1}{N} \sum_{i=1}^N (I_i - \bar{I})^2} \quad (2)$$

Where I_i is the intensity of the i -th pixel, \bar{I} is the mean intensity, and N is the total number of pixels.

- **Entropy Contrast:** Measures the randomness or complexity of the image's intensity distribution.

Formula:

$$H = - \sum_i P(i) \log_2 P(i) \quad (3)$$

Where $P(i)$ is the probability of the i -th intensity level.

2.2 Image Enhancement Techniques

- **Histogram Equalization:** Enhances the contrast of the image by spreading out the most frequent intensity values.
Result: Improved visibility of structures in the image.
- **CLAHE Enhancement:** Applies Contrast Limited Adaptive Histogram Equalization to improve local contrast.
Result: Enhanced details in specific regions of the image.

2.3 Noise Reduction Techniques

- **Gaussian Blur (Linear Filtering):** Reduces image noise by averaging pixels with a Gaussian kernel.
Formula:

$$G(x, y) = \frac{1}{2\pi\sigma^2} \exp\left(-\frac{x^2 + y^2}{2\sigma^2}\right) \quad (4)$$

Where σ is the standard deviation of the Gaussian distribution.

- **Median Filtering:** Reduces noise while preserving edges by replacing each pixel with the median of its neighborhood.
Result: Effective noise reduction with edge preservation.

2.4 Edge Enhancement Techniques

- **Laplacian Edge Enhancement:** Highlights edges by calculating the Laplacian of the image.
Formula:

$$\nabla^2 I = \frac{\partial^2 I}{\partial x^2} + \frac{\partial^2 I}{\partial y^2} \quad (5)$$

Where I is the image intensity.

- **Sobel Edge Enhancement:** Computes the gradient of the image intensity, emphasizing edges.
Formula:

$$G_x = \frac{\partial I}{\partial x}, \quad G_y = \frac{\partial I}{\partial y} \quad (6)$$

Result: Combined gradient magnitude

$$G = \sqrt{G_x^2 + G_y^2} \quad (7)$$

highlights the edges in both x and y directions.

2.5 Segmentation Techniques

- **Threshold Segmentation:** Segments the image by thresholding pixel intensities.

Formula:

$$S(i, j) = \begin{cases} 255 & \text{if } I(i, j) > T \\ 0 & \text{if } I(i, j) \leq T \end{cases} \quad (8)$$

Where T is the threshold value.

- **Watershed Segmentation:** Uses a marker-based approach to segment distinct regions in the image.

Process involves:

- Thresholding
- Noise removal via morphological operations
- Distance transformation
- Marker labelling
- Applying the watershed algorithm

3 User Interface (UI) and Interaction

3.1 Settings

- Users can upload a medical image in grayscale (PNG, JPG, JPEG).
- ROI selection allows users to focus on specific areas of the image.

3.2 Quality Measures

- Users can choose and calculate different quality measures to assess image contrast and complexity.

3.3 Image Enhancement

- Users can apply histogram equalization or CLAHE to improve image contrast.

3.4 Noise Reduction

- Gaussian blur and median filtering options are available for noise reduction.

3.5 Edge Enhancement

- Users can enhance edges using Laplacian or Sobel methods.

3.6 Segmentation

- Threshold and watershed segmentation techniques are provided for separating regions in the image.

3.7 Features

The application offers a full suite of advanced image processing functionalities to measure, enhance, and segment images. Here's what you can do:

3.7.1 1. Image Quality Measures

Assess the quality of medical images using a variety of contrast measurement techniques. This can be useful for determining image clarity and suitability for analysis.

- **Global Contrast (Michelson):** Measures the contrast based on intensity range, giving insight into overall image clarity.
- **RMS Contrast:** Computes the contrast using the Root Mean Square method, which helps to identify the spread of pixel intensity.
- **Entropy:** Calculates the complexity of pixel intensity distribution, an essential metric in assessing image information richness.

3.7.2 2. Image Enhancement Techniques

Enhancing medical images can bring out subtle details that are crucial for diagnostic accuracy. Choose from:

- **Histogram Equalization:** Redistributes pixel intensities, making hidden details more visible by balancing image brightness.
- **CLAHE (Contrast Limited Adaptive Histogram Equalization):** An advanced method that prevents excessive brightness in localized areas, ideal for bringing out fine textures in complex images.

3.7.3 3. Noise Reduction Techniques

Noise in medical images can obscure important details. The tool offers:

- **Gaussian Blur (Linear Filtering):** Smooths the image, reducing noise while maintaining essential details.
- **Median Filtering:** Reduces noise by replacing each pixel with the median of surrounding pixels, preserving edges effectively.

3.7.4 4. Edge Enhancement Techniques

Edge detection highlights structures and boundaries, aiding in medical image analysis:

- **Laplacian Edge Enhancement:** Emphasizes regions where pixel intensity changes, enhancing the edges.
- **Sobel Edge Enhancement:** Calculates gradients in both horizontal and vertical directions, allowing fine edge detection.

3.7.5 5. Segmentation Techniques

Image segmentation isolates regions of interest, making it easier to analyze specific structures:

- **Threshold Segmentation:** A quick, efficient way to separate foreground and background.
- **Watershed Segmentation:** A sophisticated technique ideal for separating overlapping structures by treating intensity as a topographic surface.

sectionStreamlit Application Layout Figure 1 shows the layout of the Streamlit app. The application interface is designed to allow users to upload images and select a contrast metric for analysis. Upon processing, the app displays the original and processed images, along with the calculated contrast metrics.

4 Example Run with Histopathological Images

Figure 2 illustrates an example run of the application using histopathological images. In this example, the uploaded image undergoes various processing steps to measure its contrast. The output includes the contrast value calculated based on the selected metric and a visualization of the processed image.

5 Explanation of the Code

The Streamlit application is implemented in Python using the Streamlit library for the user interface and OpenCV for image processing. The application workflow is as follows:

- The user uploads an image, which is displayed on the app.
- The user selects a contrast metric from a dropdown menu (e.g., Michelson contrast, RMS contrast, or Entropy).
- The app processes the image using the selected metric and displays the processed image along with the computed contrast value.

The code leverages image processing techniques to enhance contrast and provide meaningful insights into the quality of histopathological images.




Figure 1: Layout of the Streamlit application for medical image processing. The app allows users to upload histopathological images, select a contrast metric, and view the processed results.

5.1 Getting Started

Ready to dive into image processing? Here's how to set up the tool and get started.

5.1.1 Installation

1. Clone the repository:

```
git clone https://github.com/yourusername/Advanced-Medical-Image-Processing-Tool.git
cd Advanced-Medical-Image-Processing-Tool
```




Figure 2: Example output from the Streamlit application using histopathological images. The processed image demonstrates the contrast enhancement and the calculated contrast value based on the chosen metric.

2. Install dependencies:

```
pip install -r requirements.txt
```

5.1.2 Running the App

Once the setup is complete, launch the Streamlit app:

```
streamlit run MIP_FISAC.py
```

Upload a medical image in .png, .jpg, or .jpeg format and explore the options in the side panel.

5.2 How It Works

This tool provides a guided interface where you can:

1. **Upload an Image:** Start by uploading a grayscale medical image. The tool will display the image for easy reference.
2. **Select Region of Interest (ROI):** Use the sliders to define a specific region to analyze, zooming in on areas with high diagnostic value.
3. **Choose Processing Techniques:** Select from a range of quality measures, enhancement options, noise reduction, edge detection, and segmentation methods to transform the image. The tool then processes the image in real time, displaying each transformation.
4. **Download Processed Images:** Save any version of the image that's been processed for use in reports, presentations, or further analysis.

5.3 Example Workflow

Here's an example of how you might use this tool in a typical analysis workflow:

1. **Upload an MRI scan:** Open an MRI scan, adjust the region of interest to focus on a specific area.
2. **Measure Contrast:** Calculate the RMS contrast to assess the scan quality.
3. **Enhance the Image:** Apply CLAHE to bring out finer details.
4. **Reduce Noise:** Use Median Filtering to smooth out unwanted noise.
5. **Edge Detection:** Apply Sobel Edge Enhancement to highlight structural boundaries.
6. **Segmentation:** Finish with Watershed Segmentation to isolate key structures.

5.4 Technology Stack

The following libraries power this tool:

- **Streamlit:** Provides a responsive and interactive user interface.
- **OpenCV:** Performs image processing and computer vision tasks.
- **NumPy:** Enables efficient numerical operations on image data.

5.5 Future Directions

We're excited about the potential for this tool and have plans to introduce:

- **Color Image Support:** Extend functionality to color medical images, enabling analysis on more complex datasets.
- **Additional Processing Techniques:** Incorporate methods like Fourier Transform analysis and multi-scale image processing.
- **3D Visualization:** Explore support for viewing and processing volumetric (3D) medical imaging data.

5.6 License

This project is open-source and available under the MIT License.

5.7 Acknowledgments

This tool was developed with a focus on simplicity, utility, and accuracy to help the medical community leverage digital tools in image analysis. Your feedback and suggestions are always welcome!

6 Conclusion

The Advanced Medical Image Processing Tool offers a comprehensive suite of functions for analyzing and improving medical images. By integrating various processing techniques, it enhances the clarity and quality of images, aiding in better diagnosis and research. The user-friendly interface ensures easy access to powerful image processing capabilities, making it a valuable tool for medical professionals and researchers.

This detailed approach not only provides the technical foundation for each processing technique but also highlights their practical applications in enhancing medical images.